



## 24" Flammable Material Refrigerator Storage

Model # MS24RASFRW

### Standard Features

- Specially designed compressor with all electrical components fully enclosed, to prevent an electrical spark from accidentally igniting the flammable vapors from the refrigerator's contents
- Static condenser mounted to the back of the cabinet, free of electric condenser fan motor to eliminate sources of accidental ignition of flammable vapors
- Adjustable mechanical thermostat, fully enclosed to prevent ignition of flammable vapors
- Auto cycle defrost
- 2 adjustable and removable interior wire shelves
- 3 interior door shelves
- "Box Style" door
- Discrete top edge mounted door handle
- 4 leveling legs for easy height adjustment
- Jacketed commercial power cord with NEMA 5-15P plug
- Durable white powder coat paint finish standard
- Standard door lock for secure storage needs

### Compliance Specifications & Standards

- Listed in Underwriters Laboratories File SA7488(N)
- Meets National Fire Protection Association Standard 70

### Product Dimensions

23 7/8" W x 34" minimum H x 26" D (excluding door handle)



### Finish/Door Options



<b>Width</b>	23 7/8"
<b>Interior Finish</b>	White
<b>Capacity</b>	6.1 Cubic Feet; 173 Liters
<b>Toe Grill Finish</b>	Black
<b>Hinge</b>	Nylon bushing with hinge pin
<b>Handle</b>	Discrete
<b>Lock</b>	Yes
<b>Temperature Range</b>	33°F to 40°F (1°C to 5°C)
<b>Electrical Requirements</b>	115V / 60Hz / 2.5A
<b>Length of Power Cord</b>	5'
<b>Product Dimensions</b>	23 7/8" W x 34" H x 26" D
<b>Shipping Weight</b>	120 lbs; 130 lbs with skid
<b>Finishes</b>	White

### Ordering Details

Solid White Cabinet and Door with Lock and Right Hinge

MS24RASFRW